



Be Right™



Sension pH combination Gel electrode

Produktnr: 5193500
SEK Pris: Kontakta oss
Tillgänglig

Standard combination pH electrode with a non-refillable gel reference and built-in temperature sensor. Contains fixed cable with 5-pin circular connector.

Specifikationer

Diameter:	12.1 mm
Elektrolyt:	Non-refillable gel
Givar koppling:	Glass frit
Givartyp:	Glass
Kabellängd:	0.91 m
Längd:	152.5 mm
Material:	Sensor Body: Epoxy
Mätområde:	0 - 14 pH
Temperaturområde :	Continuous use: 0 - 45 °C