



CDC866T Konduktivitetscell för allmänt bruk, 4-polig, Temp.

Produktnr:	E61M015
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

4-pole conductivity electrode with MAB6 plug and epoxy body. Cell constant 1cm^{-1} . Platinised, with temperature sensor. Max. temperature 80°C , min. immersion depth 35mm.

Specifikationer

Conductivity Cell:	CDC866T
Elektrod:	4-pole, platinised
Garanti :	12 månader
Specialfunktioner:	Conductivity electrode with MAB6 plug and epoxy body.