



pH test kit, 5.5 - 8.5 pH, model 17F

Produktnr:

147006

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	200
Kapslingstyp:	B
Mätområde:	5.6 - 8.4 pH
Metodnamn:	Color disc/ Bromthymol blue
Modell:	17F
Parameter:	pH, mid range
Plattform:	Color Disc
Smallest increment:	0.1 pH
Test-kit modell:	17F