



Test kit, phosphorus, orthophosphate (reactive) , model PO-19A

Produktnr:

224801

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Simple, low-cost color disc method is the most accurate visual test kit method

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	F
Mätområde:	0 - 40 mg/L PO ₄
Metodnamn:	Color disc/Ascorbic acid
Modell:	PO-19A
Parameter:	Phosphorus
Plattform:	Color Disc
Smallest increment:	0.02
Standardmätområde:	0 - 0.8 mg/L PO ₄ ; 0 - 4 mg/L PO ₄ ; 0 - 40 mg/L PO ₄
Test-kit modell:	PO-19A

(Includes filtration for turbid samples.)