



Test kit, phosphate, model PO-23

Produktnr:

224902

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	F
Mätområde:	0.1 - 4 mg/L PO ₄
Metodnummer:	Color disc/Ascorbic acid
	for clear samples
Modell:	PO-23
Parameter:	Phosphate
Plattform:	Color Disc
Smallest increment:	0.1
	1
Standardmätområde:	0.1 - 4 mg/L PO ₄ ; 1 - 40 mg/L PO ₄
Test-kit modell:	PO-23