



Test Kit Iron, 0 - 4 mg/L Fe, Colour disc method, 100 tests

Produktnr:

146400

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	Medium Single Range
Metodnamn:	Color disc/ 1,10 Phenanthroline
Modell:	IR-18
Parameter:	Iron
Plattform:	Color Disc
Smallest increment:	0.1
Standardmätområde:	0 - 4 mg/L Fe
Test-kit modell:	IR-18
Typ:	Single Metric