



Iron color disc test kit, model IR-18A

Produktnr:

146500

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Using the simple color disc method, this kit measures Iron (Range - 0.02-1 mg/L) and contains 100 tests.

Simple, low-cost color disc method is the most accurate visual test kit method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	Low Single Range
Metodnamn:	Color disc/ 1,10 Phenanthroline
Modell:	IR-18A
Parameter:	Iron
Plattform:	Color Disc
Smallest increment:	0.02
Standardmätområde:	0.02 - 1.00 mg/L Fe
Test-kit modell:	IR-18A
Typ:	Single Metric