



## Test Kit, chlorine (free & total) & pH, model CN-67

**Produktnr:**

**1411100**

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Control of chlorine and pH in pools and spas to maintain safe bathing conditions and to protect against system corrosion. Supplied with color discs and comparator for monitoring chlorine and pH, the Free & Total Chlorine and pH Test Kit provides more accuracy than test kits using color cubes or step-wise comparators.

**One of Hach's most popular test kits**

**Includes color discs and a comparator for monitoring chlorine and pH**

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifikationer

|                |                                |
|----------------|--------------------------------|
| Antal tester : | 300                            |
| Applikation:   | Pool and Spa                   |
| Mätområde:     | 6.6 - 8.4 pH units             |
| Modell:        | CN-67                          |
| Parameter:     | Chlorine (Free & Total) and pH |
| Plattform:     | Color Disc                     |