



Hardness, iron, and pH test kit, model HA-62

Produktnr:

183700

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Model HA-62

A multi-parameter test kit for hardness, iron & pH, this kit provides an easy-to-use, cost effective way to determine the quality of your water.

Popular tests at a great price in Hach's convenient and protective test kit carrying case.

Uses color disc and drop count titration methods

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Applikation:	Water Quality
Mätområde:	6.6 - 8.4 pH units
Modell:	HA-62
Parameter:	Hardness, Iron, and pH
Plattform:	Color Disc