



Test Kit Drop Count Titrator, CA-23 Carbon Dioxide

Produktnr:

143601

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Carbon dioxide test kit with drop count titration covers ranges of 1.25-25, 2-40, and 5-100 mg/L as CO₂, with increments of 1.25, 2, and 5 respectively. Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Includes rugged carrying case, prepackaged reagent set, detailed instructions, and all apparatus you need for fast, easy, and accurate carbon dioxide test results in the field.

Specifikationer

Antal tester :	200
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	D
Mätområde:	5 - 100 mg/L CO ₂
Metodnamn:	Drop count titration/ Sodium hydroxide Phenolphthalein
Parameter:	Carbon Dioxide - As CO ₂
Plattform:	Drop Count Titration
Smallest increment:	1.25 2 5
Test-kit modell:	CA-23