



## Test kit cyanide, model CYN-3

**Produktnr:**

**201002**

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



The main sources of cyanide, an extremely poisonous substance, are metal-cleaning and electroplating baths, gas scrubber, gas works, coke ovens, and various chemical processes. It can be removed easily from waste effluents through alkaline or neutral chlorination procedures.

**Simple, low-cost method using a stepped color comparator**

**Achieve fast results with an easy-to-carry test kit**

---

### Specifikationer

|                     |                                        |
|---------------------|----------------------------------------|
| Antal tester :      | 100                                    |
| Kapslingstyp:       | R                                      |
| Mätområde:          | 0 - 0.3 mg/L CN <sup>-</sup>           |
| Metodnamn:          | Color disc/<br><br>Pyridine-Pyrazolone |
| Modell:             | CYN-3                                  |
| Parameter:          | Cyanide                                |
| Plattform:          | Color Disc                             |
| Smallest increment: | 0.02                                   |
| Standardmätområde:  | 0 - 0.3 mg/L CN <sup>-</sup>           |
| Test-kit modell:    | CYN-3                                  |