



Chloride Low Range Test Kit, Model 8-P

Produktnr:

144001

SEK Pris:
Tillgänglig

Kontakta oss

Farligt



Chloride, single parameter test kit, drop count titration (range: 5 - 100

mg/L) / silver nitrate (range: 20 - 400 mg/L), blue polypropylene case

Using silver nitrate titrant, this kit provides fast, accurate detection of chloride.

Low cost and easy to use drop count titration method using a silver nitrate titrant

Ideal for quickly determining quantitative results

Specifikationer

Antal tester :	100
Kapslingstyp:	B
Mätområde:	5 - 100 mg/L & 20 - 400 mg/L Cl ⁻
Metodnamn:	Drop count titration/ Silver nitrate
Parameter:	Chloride - As Cl ⁻
Plattform:	Drop Count Titration
Smallest increment:	5
Test-kit modell:	8-P