



Hardness (Total & Calcium) Test Kit, Model HAC-DT

Produktnr:

2063900

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifikationer

Antal tester :	100 Ca
	+100 T
Kapslingstyp:	F
Mätområde:	10 - 4,000 mg/L CaCO ₃
Metodnamn:	Digital Titrator/EDTA
Parameter:	Hardness, total & Calcium - As CaCO ₃
Plattform:	Digital Titrator
Smallest increment:	0.1-10
Test-kit modell:	HAC-DT
Vad ingår i leveransen:	Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.

Vad finns i förpackningen

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.