



Test Kit, chlorine (free and total), model CN-80

Produktnr:

2129000

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0 - 0.68 mg/L Cl ₂
Metodnamn:	Color Disc/DPD
Modell:	CN-80
Parameter:	Chlorine, low range - As free & total Cl ₂
	total Cl ₂
	(0-10 mg/L: T only)
Plattform:	Color Disc
Smallest increment:	0.02
	0.1
	0.5
Standardmätområde:	0 - 0.68 mg/L Cl ₂ ; 0 - 3.4 mg/L Cl ₂ ; 0 - 10 mg/L Cl ₂
Test-kit modell:	CN-80 ²