



Chlorine (total) test kit, model CN-21P, 10 - 200 mg/L, 100 tests

Produktnr:

2444400

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Method / Chemistry: Drop count titration / Thiosulfate

Chlorine Test Kits are manufactured by Hach Company to ensure compliance with stringent quality control standards. Using the Drop Count Titration method, the CN-21P Chlorine Test Kit is used to determine high concentrations of chlorine in sanitizing applications. 10 - 200 mg/L; contains 100 tests.

Simple, economical titration method

Kit contains all required reagents, solutions, and glassware

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	B
Mätområde:	10 - 200 mg/L
Metodnamn:	Drop count titration/ Thiosulfate
Parameter:	Chlorine, mid and high range - As total Cl ₂
Plattform:	Drop Count Titration
Smallest increment:	10
Test-kit modell:	CN-21P