



Be Right™



Dissolved oxygen test kit, model OX-2P

Produktnr:

146900

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Low dissolved oxygen levels in surface water systems usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. This test kit utilizes drop count titration and contains 100 tests.

Method / Chemistry: Drop count titration

Replacement reagent set PN 2439802

Utilizes drop count titration - a simple, economical titration method

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	0.2 - 4.0 mg/L & 1 - 20 mg/L O ₂
Metodnamn:	Drop count titration/ Modified Winkler
Parameter:	Dissolved Oxygen - As DO
Plattform:	Drop Count Titration
Smallest increment:	0.2
Test-kit modell:	OX-2P