



AccuVac[®] kit, dissolved oxygen

Produktnr:

2515050

SEK Pris:

Kontakta oss

Inget leveransdatum anges

Farligt



Low DO levels usually indicate serious pollution. Adequate amounts are crucial for fish life, but conversely, DO must be excluded from boiler feed water to prevent corrosion. The simple, low-cost color disc method is the most accurate visual test kit method.

Simple, low-cost color disc method is the most accurate visual test kit method

Achieve fast results with an easy-to-carry test kit

Specifikationer

Antal tester :	25
Kapslingstyp:	G
Mätområde:	0 - 15 mg/L O ₂
Metodnamn:	Color disc/HRDO
Modell:	Color disc/HRDO
Parameter:	Dissolved Oxygen - As DO
Plattform:	Accuvac
Smallest increment:	0.2
Test-kit modell:	AccuVac [®] Kit