



Sulfite Test Kit, Model SU-5

Produktnr:

148002

SEK Pris:
Tillgänglig

Kontakta oss

Farligt



Utilizing drop count titration/iodometric to determine sulfite, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Mätområde:	1.0 - 20 mg/L & 10 - 200 mg/L SO_3^{2-}
Metodnamn:	Drop count, titration Iodometric
Parameter:	Sulfite
Plattform:	Drop Count Titration
Smallest increment:	1 10
Test-kit modell:	SU-5