



Be Right™



Test kit, sulfate, model SF-1

Produktnr:

225100

SEK Pris:

Kontakta oss

Skickas inom 2 veckor

Farligt



Specialty Test Kits - Single Parameter

Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

Specifikationer

Antal tester :	100
Kapslingstyp:	D
Leveransvikt :	1.5
Mätområde:	50 - 200 mg/L
Metodnamn:	Extinction/Turbidimetric
Parameter:	Sulfate - As SO_4^{2-}
Smallest increment:	50
Test-kit modell:	SF-1