



Test kit, alkalinity, model AL-AP

Produktnr:

2444300

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Utilizing drop count titration/sulfuric acid to determine alkalinity, this kit contains 100 tests.

Utilizes drop count titration - a simple, economical titration method

Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm. **1 gpg = 17.1 mg/L (MO) = Methyl Orange
Kapslingstyp:	B
Mätområde:	0 - 8 gpg & 1 - 20 gpg
Metodnamn:	Drop count titration/Sulfuric acid
Parameter:	Alkalinity - (P) & Total (MO) As CaCO ₃
Plattform:	Drop Count Titration
Smallest increment:	0.4 gpg
	1 gpg
Test-kit modell:	AL-AP