



Alkalinity test kit AL-AP, 5 - 100 / 20 - 400 mg/L CaCO₃ 100 pcs

Produktnr:

2444301

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists.

Method: Drop count titration

Caused by the presence of carbonates, bicarbonates, hydroxides, and other dissolved salts, alkalinity is important to food and beverage processors, boiler operators, and aquatic biologists. Kit utilizes drop count titration method and contains 100 tests.

Simple, economical titration method

Rugged kit contains all required reagents, solutions, and glassware

Specifikationer

Antal tester :	100
Fotnot :	*mg/L unless otherwise noted, ppb = µg/L, ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Kapslingstyp:	B
Mätområde:	5 - 100 mg/L & 20 - 400 mg/L
Metodnamn:	Drop count titration/Sulfuric acid
Parameter:	Alkalinity - (P) & Total (MO) As CaCO ₃
Plattform:	Drop Count Titration
Smallest increment:	5
	20
Test-kit modell:	AL-AP, MG-L