



Bottle, BOD, 300 mL, 24/pk

Produktnr:	62124
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

The standard for BOD analysis. Borosilicate clear glass bottles with flared lip and ground-glass stoppers to prevent air entrapment and for water seal. Frosted marking area. Stoppers included. Pack of 24.

Specifikationer

Body Material:	Glass
Kapacitet :	300 mL
Material:	Borosilicate Glass
Package size UOM:	24pk
Volym/förpackningsstorlek:	24 /pk