



Beaker, PP, low form, 100 mL

Produktnr:	108042
SEK Pris:	Kontakta oss
Tillgänglig	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Bull nose lip and molded pour spout minimise spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikationer

Diameter, yttre:	63.0 mm with spout
Enhet:	1 each
Kapacitet :	100 mL
Material:	Polypropylene, Low Form
Volym/förpackningsstorlek:	1 each