



pH electrode 5338, plastic, gel, fixed cable, 0-80 °C, 2 bar

Produktnr:	LZU5338.97.0002
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Low cost pH electrode 5338 with plastic body, gel electrolyte and fixed cable. 0 - 14 pH, 0 - 80 °C, up to 2 bar.

Specifikationer

Elektrolyt:	Gel
Garanti :	6 månader
Material:	Plastic
Mätområde pH:	0 - 14 pH
Parameter:	pH
Temperaturområde :	0 - 80 °C
Tryckområde:	Up to 2 bar