



5334T, pH-elektrod med Pt1000, 0–100 °C, 5 bar, MP-4-kontakt

Produktnr: LZU5334T.97.002

SEK Pris: Kontakta oss

Tillgänglig

pH electrode with integrated PT1000 temperature sensor, gel electrolyte and MP-4 plug. 0 - 14 pH, 0 - 100 °C, up to 5 bar. Ring of porous PTFE.

Applications: waste water,

galvanising industry, osmosis processes, demineralised water, dirty or viscous samples, etc.

Specifikationer

Diameter:	12 mm
Elektrolyt:	Gel
Garanti :	6 månader
Kabelanslutning :	MP-4
Längd:	120 mm
Material:	Givare: glas
Mätområde pH:	0 - 14 pH
Membran:	PTFE-ring
Parameter:	pH
Reference system:	Inkapslade Ag/AgCl-kristaller
Temperaturgivare:	Pt1000
Temperaturområde :	0 - 100 °C
Tryckområde:	upp till 5 bar