



pH electrode 5364, for difficult media, 0-130 °C, 6 bar, sterilisable, S8 screw head

Produktnr:	LZU5364.97.0002
SEK Pris:	Kontakta oss
Tillgänglig	

Sterilisable pH electrode 5364 with polymer electrolyte and S8 screw head. 0 - 14 pH, 0 - 130 °C, up to 6 bar. Open diaphragm.

Specifikationer

Diameter:	12 mm
Elektrolyt:	Polymer
Garanti :	6 månader
Kabelanslutning :	S8
Längd:	120 mm
Material:	Sensor: glass
Mätområde pH:	0 - 14 pH
Membran:	Open, 2 laterals
Parameter:	pH
Reference system:	Encapsulated Ag/AgCl crystals
Temperaturområde :	0 - 130 °C
Tryckområde:	upp till 6 bar