



## Conductivity cell 5390 with Pt1000, k=0.1, 100°C, MP-5, 1/2"

|                           |                        |
|---------------------------|------------------------|
| <b>Produktnr:</b>         | <b>LZU5390.97.0002</b> |
| SEK Pris:                 | Kontakta oss           |
| Inget leveransdatum anges |                        |

Conductivity cell 5390 + PT1000 A.T.C., for pipe, in titanium, PVDF and stainless steel, 0,1...10000 µS/cm, 0...100 °C. MP-5 connector, ½" thread.

---

### Specifikationer

|                         |                            |
|-------------------------|----------------------------|
| Cellkonstant k:         | 0.1 cm <sup>-1</sup>       |
| Kabelanslutning :       | MP-5                       |
| Längd:                  | 50 mm                      |
| Material:               | Measuring Cell: Titanium   |
| Mätområde:              | 0.10 - 10.000 mg/L (µS/cm) |
| Temperaturkompensation: | 0 - 100 °C                 |
| Temperaturområde :      | Pt 1000                    |