



Be Right™



EVITA OXY 1100 Membrane 125µm for Dissolved Oxygen sensor, 1 piece

Produktnr: 085G0023
SEK Pris: Kontakta oss
Tillgänglig

Recommended for measuring range 2-50mg/L DO.

Specifikationer

Avfall:	OXY 1100 sensor must be disposed as ordinary industrial waste.
Mätområde:	2.0 - 50 mg/L (ppm) recommended value to maintain operating life
Responstid :	110 s
Specialkommentarer:	This sensor is a part of EVITA OXY 4100 and OXY 4150 sensor.
Tjocklek, membran:	125 µm