



Digital Titrator-patron med silvernitrat, 0.2256 N

Produktnr:

1439601

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Bruksfärdig titreringslösning för bestämning av klorid (silvernitratmetoden), 10-10000 mg/l Cl⁻.

For chloride determination by Digital Titrator argentometric titration. Hach method 8207. Digits x 0.1 = mg/L Cl⁻ in 100 mL sample.
Approximately 40-100 tests

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
Mätområde:	10 - 100 mg/L Cl ⁻
Metodnamn:	Silver Nitrate
Metodnummer:	8207
Parameter:	Chloride
Plattform:	Digital Titrator
Titrant:	0.2256 N AgNO ₃