



Natriumhydroxid Digital Titrator (patron), 1.600 N

Produktnr:

1437901

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Bruksfärdig titreringslösning för syra-bastitrering (1-10 meq/ syra) och bestämning av aciditet (100-4000 mg/l CaCO₃)

For acidity determination by Digital Titration, 100 - 4000 mg/L CaCO₃

Accurate results - when and where you need them

Digital Titrator methods use a high-precision dispensation device and cartridges of concentrated titrant to allow analyses in the lab, the plant or the field.

Fast and Flexible

Interchangeable cartridges allow you to change titrants in seconds. No cleaning or rinsing burets.

Versatile

More than 50 parameter/range combinations are available

Economical

Concentrated titrants allow 40-100 typical titrations per cartridge.

Specifikationer

Antal tester :	40 - 100
EPA compliant:	Ja
Mätområde:	1 - 4000 mg/L CaCO ₃
Metodnamn:	Fenolftaleína
Metodnummer:	8201 (MO)
Parameter:	Acidity
Plattform:	Digital Titrator
Titrant:	1.600 N NaOH