



ULR-verifieringssats, Stablcal, 100 mL

Produktnr:

2714600

Farligt

SEK Pris:

Kontakta oss

Inget leveransdatum anges



Standardpaket för kontroll av turbiditet i låga mätområden, inkluderar 100 ml av följande: 0,10 NTU, 0,30 NTU, 0,50 NTU, 1,0 NTU standard, silikonolja och oljedukar samt instruktioner.

Stablcal Verification Kit includes: 0.10 NTU, 0.30 NTU, 0.50 NTU, 1.0 NTU standards (100 mL each), silicone oil and oiling cloth, and instructions. For verification only. Not for calibration.

Precise, Pre-mixed Formazin Solutions

With proprietary manufacturing technology, Hach prepares Stablcal Standards in precise concentrations at a wide variety of concentrations - providing accurate turbidimeter calibrations.

Based on Formazin Primary Standards

Stablcal Stabilized Formazin Standards contain the same light scattering polymer as traditional Formazin primary turbidity standards - and perform equivalently.

Enhanced Stability Saves Time and Labor

While Formazin primary standards can be prepared directly in the lab, the process is labor intensive and time consuming - and the diluted standards are not stable. Stablcal is a stabilized version of Formazin - providing repeatable results, without the time and effort required for lab-prepared standards.

Confidence in results

Hach was the first to apply Formazin as a turbidity standard in the drinking water industry, more than 40 years ago. Stablcal standards are based on decades of experience in turbidity measurement.

Specifikationer

Fotnot :	*For verification only. Not for calibration.
Koncentration:	0,1, 0,3, 0,5, 1,0 NTU
Volym/förpackningsstorlek:	paket med 4