



SP510 Hardness Buffer Solution (0.3, 1, 2, 5 mg/L) Hardness buffer solution, 500mL

Produktnr:

2768549

Farligt

SEK Pris:

Kontakta oss

Skickas inom 2 veckor



For use with the SP510 Hardness analyser.

Specifikationer

Instrument:	SP510
Koncentration:	0.3, 1, 2, 5 mg/L
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Parameter:	Hardness