



Amino Acid F Reagent for Silica analysis, 2.9 L

Produktnr:

2353103

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Amino Acid F Reagent (2.7 L) used in the Series 5000 Silica Analyser.

Amino Acid F reagent is supplied in two parts to ensure a long shelf life: a 3-liter bottle of Amino Acid F Solution and a bottle of Amino Acid F Powder. Powder must be added to the solution immediately before installation in the analyser.

Specifikationer

Instrument:	5500sc Silica
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Mätområde:	N/A
Parameter:	Silica
Volym/förpackningsstorlek:	1 x 2.7 L & 1 x 308 g