



Stablcal Standardlösung för turbiditet 20 NTU (1000 mL)

Produktnr:

2660153

Farligt

SEK Pris:

Kontakta oss

Tillgänglig



Stablcal Turbidity Standard for all Turbidimeters. USEPA accepted for calibration purposes. Stablcal Stabilised Formazin Turbidity Standards are true Formazin dilutions developed for use in any turbidimeter. NTU Value: 20.0. Bottle: 1000 mL (1 L).

Specifikationer

Koncentration:	20 NTU
Lagringsförhållanden:	Store between 5 - 25 °C, away from direct sunlight
Mätområde:	N/A
Parameter:	Turbidity
Volym/förpackningsstorlek:	1 L