



CDC511T Konduktivitetscell, 4-polig, rent vatten

Produktnr:	E61M009
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

4-pole conductivity electrode with MAB6 plug and PMP body. Cell constant 1cm^{-1} . With temperature sensor. Max. temperature 80°C, min. sample volume 3mL. For pure water.

Specifikationer

Conductivity Cell:	CDC511T
Elektrod:	4-pole, pure water
Garanti :	12 månader
Specialfunktioner:	Conductivity electrode with MAB6 plug and PMP body.