



## **CDC861T Konduktivitetsscell, temp., CDC861T Konduktivitetsscell, temp,**

|                           |                |
|---------------------------|----------------|
| <b>Produktnr:</b>         | <b>E61M016</b> |
| SEK Pris:                 | Kontakta oss   |
| Inget leveransdatum anges |                |

For use in strong acids and bases, salinity.

### **The right conductivity cell for your application.**

It is important to choose a conductivity cell with the right construction and geometry for your particular application and working conditions. The Radiometer Analytical line offers a variety of conductivity cells to meet a large number of applications.

---

### **Specifikationer**

|                    |                                                 |
|--------------------|-------------------------------------------------|
| Conductivity Cell: | CDC861T                                         |
| Elektrod:          | 4-pole, platinised                              |
| Garanti :          | 12 månader                                      |
| Specialfunktioner: | Conductivity cel with MAB6 plug and glass body. |