



pHC2011-8 Kombinerad pH-elektrod, alkal. Prover, röd stav, BNC

Produktnr: E16M317
SEK Pris: Kontakta oss
Inget leveransdatum anges

Long-term accuracy and reproducibility with Red Rod technology

pH electrode for alkaline samples. With Red Rod reference system. Measuring range 0-14 pH, 0-100 °C. Liquid junction: porous pin. With BNC plug.

The symmetry of Red Rod electrodes means that the iso pH is the same as the zero pH which provides highly reproducible results even when temperature fluctuations occur.

Reproducible measuring conditions ensured by a special bayonet head

Wide temperature range from -10 to 100°C

Fast and accurate response

Stable signal

Trouble-free liquid junction

Specifikationer

Elektrodtyp:	pHC2011-8
Elektrolyt:	Refillable: Saturated KCl
Garanti :	12 månader
Givar koppling:	Porous pin
Givartyp:	Laboratory combination pH electrode with Red-Rod reference system
Kabellängd:	1 m
Längd:	103 mm
Material:	Glass
Mätområde:	0 - 14 pH
Specialfunktioner:	Special glass membrane for alkaline samples and Red-Rod reference system
Temperaturområde :	0 - 100 °C
Termistor:	None (external required)

Obligatoriska tillbehör

- Påfyllningslösning, referens, mättad KCl (KCl.L), 100 mL (Item S21M002)
- Påfyllningslösning, referens, KCl-kristaller (KCl.C), 15 g (Item S21M001)