



XR300 AG/AGCL REF ELEKTROD

Produktnr:	B20B300
SEK Pris:	Kontakta oss
Tillgänglig	

For application of Voltammetry.

Specifikationer

Dimension (D x L):	8 mm x 120 mm
Elektrolyt:	Saturated KCl + saturated AgCl
Garanti :	6 månader
Givar koppling:	Porous pin
Givartyp:	Laboratory reference electrode
Längd:	120 mm
Material:	Glass
Temperaturområde :	-5 - 80 °C