



## Tensette Pipette 1.0 - 10.0 mL

<b>Produktnr:</b>	<b>1970010</b>
SEK Pris:	Kontakta oss
Tillgänglig	

Perform routine laboratory pipetting quickly, economically, and accurately with Hach's automatic Tensette Pipette. Two models dispense from 0.1 - 1 mL or from 1 - 10 mL in 0.1 or 1.0 mL increments, respectively. The pipette saves time, reduces glassware cleanup and assures accuracy. A Set includes one pipette, a manual and 50 disposable tips.

---

### Specifikationer

Volym/förpackningsstorlek:	1 each
----------------------------	--------