



TMS Filtration Probe, complete incl. probe body

Produktnr:	LZH160
SEK Pris:	Kontakta oss
Tillgänglig	

Specifikationer

Diameter:	45 mm
Längd:	485 mm
Material:	Membran housing: PVC
	Gaskets: FKM/FPM
	Membran: Porospoly (pore size 2 µm)
Omgivningstemperatur:	4 - 40 °C
Provvolyt:	200 - 500 mL/h, depending on pollution level of the membrane
Vikt :	0,9 kg