



pHRETPT100 pH sensor, 12 mm, retractable

Produktnr:	LZX477
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

With integrated temperature sensor PT100 and 10 m cable. Suitable for pressure lines.

Specifikationer

Elektrolyt:	KCl/AgCl + KNO ₃ gel, suitable for high temperatures
Garanti :	6 månader
Kabellängd:	10 m fixed cable
Längd:	225 mm
Material:	Stainless steel
Mätområde:	0 - 14 pH
Membran:	Porous PTFE, double ring diaphragm
Neddopningsdjup:	205 mm
Omgivningstemperatur:	-5 - 135 °C
Reference system:	Ag/AgCl
Temperaturgivare:	PT100
Tryckområde:	15 bar