



Replacement pH electrode for digital 1200-S sc pH sensor

Produktnr:

LZX889

SEK Pris:

Kontakta oss

Skickas inom 3 dagar

Specifikationer

Diameter:	12 mm
Garanti :	6 månader
Kabelanslutning :	S7 plub connector
Känslighet:	± 0,5 %
Material:	Sensor: Glass (general purpose)
Noggrannhet:	± 0.02 pH
Repeterbarhet:	± 0.5 %
Responstid T90:	< 15 s
Specialkommentarer:	Only applicable in combination with 1200-S sc pH sensor