



AD digital gateway for 1200 sc with cable gland kit

Produktnr:

LZY328

SEK Pris:

Kontakta oss

Tillgänglig

Gateway to connect analogue 1200 pH and ORP sensors to SC controllers

Specifikationer

Diameter:	34 mm
Elektrodkompatibilitet:	All combination pH and ORP electrodes with and without solution ground
Garanti :	24 månader
Kommunikation:	MODBUS
Längd:	175 mm
Omgivningstemperatur:	-20 - 60 °C at 0 - 95% RH (non-condensing)
Specialkommentarer:	Cable gland kit composed of: Strain relief M20 x 1.5 seal inserts for different cable diameters instruction sheet DOC273.99.00052
Vikt :	0,145 kg