



Pole mounting hardware pH, 10cm bracket, SS pole 2m

Produktnr: LZY714.99.23920
SEK Pris: Kontakta oss
Tillgänglig

Stainless steel pole mounting hardware for 1" NPT pH sensor including 10 cm SS bracket and 2 m SS pole. Suitable for all 1" NPT pH sensors like 1200S sc (LXV426.99.10001) and pHD ORP (LXV427.99.10001)

Specifikationer

Diameter:	40 mm
Instrument:	suitable for: all 1"NPT pH sensors.
Längd:	Pole: 2000 mm
Material:	Wetted parts: pole SS 1.4301; sensor adapter PVC; screw set SS.
Vikt :	Lifting part only: 4,9 kg