



Maintenance Kit for hydrogen thermal conductivity sensors with protection cap 29142 for use with the 3654

Produktnr:	32747
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Includes membranes 2952A, membrane holding ring 29228.01, O-rings 29039.0 and 28613.0, and tools for sensor maintenance.

Specifikationer

Instrument:	Orbisphere 3654 Portable N2 and H2 Analyzer (TC), Orbisphere 3658 Analyzer Portable CO2 Analyzer (TC)
Lagringsförhållanden:	Store in a dry place at room temperature
Parameter:	H2