



Magnesium Sulfate Standard, 29 mL, for Calcium Titration and Hardness Monitor, 10,000 mg/L as CaCO₃, for SP510

Produktnr:	102233
SEK Pris:	Kontakta oss
Tillgänglig	

For checking Hardness Monitor response: 10,000 ±1000 mg/L as CaCO₃

Specifikationer

Instrument:	SP510
Koncentration:	10,000 mg/L
Lagringsförhållanden:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.
Parameter:	Hardness

Obligatoriska tillbehör

- Skyddshandskar storlek L, blå, nitril, puderfria, 50 styck (Item SM743L)