



Beaker, Polypropylene, Low Form, 250 mL, 6/pk

Produktnr:	108076
SEK Pris:	Kontakta oss
Skickas inom 2 veckor	

Griffin style (low form) beakers from Nalgene are the industry standard for a wide variety of lab applications. Graduations are high quality screen printing in sizes ≥ 250 mL. Bull nose lip and molded pour spout minimize spillage in liquid transfer. Autoclavable, non-toxic and biologically inert. Volumes and graduations are approximate.

Specifikationer

Diameter, yttre:	83.6 mm with spout
Enhet:	6 /pk
Kapacitet :	250 mL
Material:	Polypropylene, Low Form