



Flask, volumetric, polypropylene, 250 mL

Produktnr: 1406046

UTGÅNGEN ARTIKEL

Den här artikeln är inte längre tillgänglig.

Polypropylene flask and screw cap. Durable with excellent chemical resistance, with the same accuracy as glass. Permanent color graduation line with marking spot. Calibrated “to contain/deliver” at 20°C. The 200/203 capacity flask has dual marks useful in DO determination via Winkler Titration.

Specifikationer

Certifieringar :	Not Applicable
Enhet:	1 each
Kapacitet :	250 mL
Material:	Polypropylene