



Tube, culture, 16 x 150 mm, screw cap, 10/pk

Produktnr:	1436120
SEK Pris:	Kontakta oss
Inget leveransdatum anges	

Cap has a firmly bonded rubber liner providing a sure and highly inert sealing surface.

Specifikationer

Diameter:	16 mm
Dimension (D x L):	16 mm x 150 mm
Kapacitet :	Capacity to Neck 20 mL
Längd:	150 mm
Specialfunktioner:	Tube, 16 x 150 mm
Volym/förpackningsstorlek:	10 /pk