



Indikator för pulverkuddar, bromphenolblått, pk/100

Produktnr: 1455099
SEK Pris: Kontakta oss
Skickas inom 2 veckor

pH indicator for acidity and alkalinity titrations. Visual transition interval: pH 3.0 (yellow)-pH 4.6 (blue). Sample Size: 50 mL

Specifikationer

Antal tester :	100
Mätområde:	0 - 5000 CaCO ₃ (MO)
Metodnamn:	Titration
Metodnummer:	8219 (MO)
Plattform:	Powder Pillows
Volym:	100 Powder Pillows